Fingertip Oximeter User Manual

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Date: _____

QA: _____

This product has been inspected in accordance with
the standards specified in the User Manual.

Shenzhen Creative Industry Co., Ltd

Instructions to User

Dear Users,

Thank you very much for purchasing our product. Please read the manual very carefully before using this device. Failure to follow these instructions can cause measuring abnormality or damage to the Oximeter.

The manual is published in English and we have the ultimate right to explain the Manual. No part of this manual may be photocopied, reproduced or translated into another language without the prior written consent. We reserve the right to improve and amend it at any time without prior notice.

Version of the Manual: Ver 1.3

Revised date: July 12, 2019

All rights reserved.

Notes

- The contents contained in this manual are

subject to change without notice.

- Information furnished by our company is believed to be accurate and reliable. However, no responsibility is assumed by us for its use, or any infringements of patents or other rights of third parties that may result from its use.

3502-2290064

Instructions for Safe Operations

- ☛ Check the device to make sure that there is no visible damage that may affect user's safety or measurement performance with regard to sensors and clips. It is recommended that the device should be inspected minimally once a week. When there is obvious damage, stop using the device.
- ☛ Special attention should be paid while the oximeter is used constantly under the ambient temperature over 37°C, burning hurt may occur because of over-heating of the sensor at this situation.
- ☛ Necessary maintenance must be performed only by qualified service technicians. Users are not permitted to maintain it by themselves.
- ☛ The Oximeter cannot be used together with devices not specified in User Manual.

Cautions

- 💣 Explosive hazard—**DO NOT** use the Oximeter in environment with inflammable gas such as some ignitable anesthetic agents.
- 💣 **DO NOT** use the Oximeter while the testee is under MRI or CT scanning.

Warnings

- 💣 An uncomfortable or painful feeling may appear if using the Oximeter continuously on the same place for a long time, especially for poor microcirculation patients. It is recommended that the Oximeter should not be applied to the same location for longer than 2 hours. If any abnormal condition is found, please change the position of Oximeter.
- 💣 **DO NOT** clip this device on edema or tender tissue.

-
- ☛ The light (the infrared light is invisible) emitted from the device is harmful to the eyes, so service technician or testee should not stare at the light.
 - ☛ The oximeter is not a treatment device.
 - ☛ The local law must be followed when disposing of the device.

Attentions

- 🔔 Keep the Oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture.
- 🔔 The device should be kept out of the reach of children.
- 🔔 If the Oximeter gets wet, please stop using it and do not resume operation until it is dry. When it is carried from a cold environment to a warm and humid environment, please do not use it immediately.
- 🔔 **DO NOT** operate the button on the front

panel with sharp materials.

 **DO NOT** use high temperature or high pressure steam disinfection on the Oximeter. Refer to Chapter 9 for instructions of cleaning and disinfection.

Declaration of Conformity

The manufacturer hereby declares that this device complies with the following standards:

IEC 60601-1

ISO 9919

and follows the provisions of the council directive MDD93/42/EEC.

Caution: U.S. federal law restricts this device to sale or use by or on the order of a physician.

1 Overview

1.1 Appearance


Figure 1 Front View

Figure 2 Rear View

1.2 Name and Model

Name: Fingertip Oximeter

Model: POD-1/POD-2/POD-3/POD-4

PC-60B1/PC-60B2/PC-60B3/PC-60B5/

PC-60C1/PC-60C2/ PC-60D2/PC-60N

1.3 Intended Use

This Fingertip Oximeter is intended for measuring the pulse rate and functional oxygen saturation (SpO_2) through patient's finger. It is applicable for spot-checking SpO_2 and pulse rate of adult and pediatric patients in homes and clinics.

2 Battery Installation


Figure 3 Battery Installation

1. Refer to Figure 3, insert two AAA size batteries into the battery compartment properly.
 2. Replace the cover.
- Please make sure that the batteries are correctly installed, or incorrect installation may cause the device not to work.

3 Quick Guide to Operation

Step	Operation for POD-1/POD-2/POD-3/POD-4, PC-60B1, PC-60B2, PC-60B3, PC-60B5, PC-60C2 and PC-60D2.
1	Open the clip and put finger into the rubber cushions of the clip (make sure the finger is in the correct position), and then clip the finger.
2	Wait 2 seconds, the Oximeter will power on automatically and start to measure.
-	Short time press Display Key to change display direction.
-	Longtime press Display Key to shift the parameter display between PR and PI. Note: POD-1/2/3 do not have PI display function.
3	Remove your finger, the Oximeter will power off automatically.

Step	Operation for PC-60C1, PC-60N
1	Open the clip and put finger into the rubber cushions of the clip (make sure the finger is in the correct position), and then clip the finger.
2	Wait 2 seconds, the Oximeter will power on automatically and start to measure.
-	Tilt the Oximeter to change display direction.
-	Short time press Display Key to shift the parameter display between PR and PI.
-	Longtime press Display Key to enter setup menu screen.
-	On setup menu screen: Short time press Display Key to choose the setting item; Longtime press Display Key to active the setting item, then short time press it to modify the setting parameter; Next, longtime press Display Key to confirm the modification and exit from this setting item.

	At last, move the setting item to “Save, exit menu”, and long time pressing Display Key to store the modification and exit from the setup menu.
3	Remove your finger, the Oximeter will power off automatically.

4 Display Screen

For PC-60B5 Oximeter, SpO₂ and PR symbols are imprinted on the front panel and its corresponding value displays on display screen below the symbol. (see figure 5).

Short time press display key, display screen (including values, pulse intensity histogram, pulse beat symbol and waveform) will be flipped 180°. (See figure 6).


Figure 5


Figure 6

Note: the two figures both show that SpO₂ value is 99% and PR value is 86 bpm.


For other type Oximeters


Measuring screen
(POD-1-2-3)

Default measuring screen
(POD-4, PC-60B1/B2/B3/B5/C1/C2/D2/N)

SpO₂&PI display screen
(POD-4, PC-60B1/
B2/B3/B5/C1/C2/D2/N)

Setup menu screen
(PC-60 C1/N)

Note: the "Setup menu screen" of PC-60C1/N shows the factory default limits setting value which is different from other models.

5 Instructions for Operation

- The finger should be put in properly and correctly.
- Do not shake the finger. Keep at ease during measurement.
- Do not put wet finger directly into sensor.
- Avoid placing the device on the same limb which is wrapped with a cuff for blood pressure measurement or during venous infusion.
- Do not let anything block the emitting light from device.
- Vigorous exercise and electrosurgical device interference may affect the measuring accuracy.
- Nail polish may affect the measuring accuracy, and too long fingernail may cause failure of measurement or inaccurate result.

-
- If the first reading appears with poor waveform (irregular or not smooth), then the reading is unlikely true, the more stable value is expected by waiting for a while, or a restart is needed when necessary.

Note: Due to the working principle of orientation sensor used in Oximeter, such as the model PC-60B3, there is a small metal ball which is movable within its compartment of the orientation-sensor. Therefore you can hear a slight “clatter” sound when you wave or shake the oximeter. It is normal and not caused by unwanted part.

6 Technical Specifications

A. SpO₂ Measurement

Transducer: dual-wavelength LED sensor
with wavelength:

Red light: 663 nm, Infrared light: 890 nm.

Maximal average optical output power:
≤1.5mW

SpO₂ measuring range: 35%~100%

SpO₂ measuring accuracy:

$A_{rms} \leq 3\%$ for SpO₂ range from 70% to 100%

Note: A_{rms} is defined as root-mean-square value of deviation according to ISO 9919.

B. Pulse Rate measurement

PR measuring range: 30bpm~240bpm

PR measuring accuracy: $\pm 2\text{bpm}$ or $\pm 2\%$
(whichever is greater)

C. Perfusion Index(PI) Display range

0%~20%

D. Preset over-limits

SpO₂ low limit: 90%

Pulse Rate: high limit: 120bpm
low limit: 50bpm

**E. Over-limit indication settings for
PC-60C1/N:**

SpO₂:

low limit setting range: 85%~95%

Default setting: 85%

Pulse Rate:

Low limit setting range: 30~60bpm;
High limit setting range: 100~240bpm;
Default setting: high: 120bpm; low: 50bpm

F. Audible & visual alert function

When measuring, if SpO₂ value or pulse rate value exceeds the preset limit, the device will alert automatically and the value which exceeds limit will flash on the screen.

G. Power supply requirement:

2 x LR03 (AAA) alkaline batteries
Supply voltage: 3.0VDC
Operating current: $\leq 40\text{mA}$

H. Environment requirement

Operating Temperature: 5°C ~40°C
Operating Humidity: 30%~80%
Atmospheric pressure: 70kPa~106kPa

I. The performance under low perfusion condition

The accuracy of SpO₂ and PR measurement still meet the precision described above when

the modulation amplitude is as low as 0.6%.

J. Resistance to interference of surrounding light:

The difference between the SpO₂ value measured in the condition of indoor natural light and that of darkroom is less than $\pm 1\%$.

K. Dimensions: 66 mm (L) $\times$ 36 mm (W) $\times$ 33 mm (H)

Net Weight: 60g (including batteries)

L. Classification

The type of protection against electric shock: Internally powered equipment.

The degree of protection against electric shock: Type BF applied parts.

Electro-Magnetic Compatibility: Group I, Class B

7 Accessories

- A. A lanyard
- B. Two batteries
- C. A pouch (optional)
- D. An External SpO₂ Probe (optional)
- E. A User Manual
- F. Quality Certificate

Note: The accessories are subject to change. Detailed items and quantity see the Packing List.

8 Repair and Maintenance

8.1 Maintenance

The service life (not a warranty) of this device is 5 years. In order to ensure its long service life, please pay attention to the maintenance.

- Please change the batteries when the low-voltage indicator lightens.
- Please clean the surface of the device before using. Wipe the device with alcohol first, and then let it air dry or wipe it dry.
- Please take out the batteries if the oximeter

will not be used for a long time.

- The recommended storage environment of the device:
ambient temperature: $-20^{\circ}\text{C} \sim 60^{\circ}\text{C}$, relative humidity 10%~95%, atmospheric pressure: 50kPa~107.4kPa.
 - The oximeter is calibrated in the factory before sale, there is no need to calibrate it during its life cycle. However, if it is necessary to verify its accuracy routinely, the user can do the verification by means of SpO_2 simulator, or it can be done by the local third party test house.
- ⚠ **High-pressure sterilization cannot be used on the device.**
- ⚠ **Do not immerse the device in liquid.**
- ⚠ **It is recommended that the device should be kept in a dry environment. Humidity may reduce the life of the device, or even damage it.**

8.2 Cleaning and Disinfecting Instruction

- Surface-clean sensor with a soft cloth by wetting with a solution such as 75% isopropyl alcohol, if low-level disinfection is required, use a 1:10 bleach solution.
- Then surface-clean with a cloth saturated with clean water and dry with a clean, soft cloth.

Caution: Do not sterilize by irradiation steam, or ethylene oxide.

Do not use the sensor if it is damaged.

9 Troubleshooting

Trouble	Possible	Solution
Display direction doesn't change or changes insensitively.	Maybe the oximeter is not used for a long time, the movable metal ball within the orientation-sensor can not move freely.	Please shake the oximeter with a certain force to make the movable metal ball move freely. If the problem still exists, maybe the orientation-sensor is not working properly. Please contact the local service center.

Trouble	Possible Reason	Solution
The SpO₂ and Pulse Rate display instable	1. The finger is not placed far enough inside. 2. The finger is shaking or the patient is moving.	1. Place the finger correctly inside and try again. 2. Let the patient keep calm.
Can not turn on the device	1. The batteries are drained or almost drained. 2. The batteries are not inserted properly. 3. The device is malfunctioning.	1. Change batteries. 2. Reinstall batteries. 3. Please contact the local service center.

**No
display**

1. The device will power off automatically when it gets no signal for 8 seconds.
2. The batteries are almost drained.

1. Normal.
2. Change batteries.

Appendix

I Key of Symbols

Symbol	Description
%SpO ₂	Pulse oxygen saturation
 BPM/PR	Pulse rate (beats per minute)
PI%	Perfusion Index (%)
	Pulse intensity bargraph
 /	Low battery voltage

CE mark

SN

Serial number

Date of manufacture

Authorised representative in the European community

Manufacturer
(including address)

With Type BF applied part

Warning — See User Manual

Disposal of this device
according to WEEE regulations

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