FR

NEOJELLY®
Gel lubrifiant stérile

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Jiangsu Province JianErKang
Medical Dressing Co., Ltd.
NO. 1, Jianerkang Road
Zhixi Industry Concentration Area, Zhixi Town
213251 Jintan Changzhou City, China

Importé par :

ASEPT INMED
9, Avenue Mercure 31130 Quint Fonsegrives France
Tél : +33(0)5 62 57 69 00 Fax : +33(0)5 62 57 69 01
www.aseptinmed.fr

Shanghai International Holding Corp. GmbH
Eiffestrasse 80, 20537 Hamburg, Germany

EN

NEOJELLY®
Sterile lubricating jelly

Please read this leaflet carefully before using this medical device because it contains important information.
If you have more questions, please contact Asept InMed.

1. DESCRIPTION

NeoJelly® devices are sterile lubricating jelly. They are available in pouch and tube.
This device is a clear, greaseless and water-soluble sterile lubricating jelly.

2. INDICATIONS

Intended use for medical purpose to lubricate body orifices to facilitate the entry of diagnostic and therapeutic devices like **insertion catheter, endoscopy, gynecological examination, rectal procedure, cystoscopy.**

3. CONTRAINDICATIONS

No known contraindication.

4. WARNINGS

- The product should be used only by qualified healthcare professionals.
- The manufacture will not be responsible for any direct incidental or consequential damages resulting from use of product not indicated in this IFU.

5. USE PRECAUTIONS

- Do not use if package is opened or damaged.
- Check the expiry date before use.
- Do not re-sterilize.
- If irritation occurs, discontinue use and consult a doctor.
- Avoid eye and ear contact.
- Keep out of the reach of children.
- Should not be used on opened wounds.

6. HYGIENE PRECAUTIONS

Respect the general rules of hygiene and asepsis of health establishments.
It is recommended to wash your hands before and after use.

7. FREQUENCY OF USE

Single use for pouches

8. MODALITY OF USE

It is recommended to wash or disinfect your hands before and after use.

Pouch:

1. Remove pouch from box;
2. Tear pouch at one end;
3. Squeeze lubricating jelly out as needed and apply.

Tube :

1. Remove cap and peel foil from end of tube;
2. Squeeze lubricating jelly out as needed and apply to instrument or skin;
3. Close after use.

9. STORAGE AND TRANSPORT CONDITIONS

These devices must be kept, stored and handled, avoiding excessive heat, protecting from direct sunlight and moisture.
Storage between 10 to 35°C.

10. DISPOSABLE CONDITIONS

Used devices must be eliminated according to the hygiene protocol of the health facility and to the actual regulation.

11. PERFORMANCES

Not applicable.

12. MATERIAL

These devices are a mixture of glycerin, polyethylene glycol carbomer, methyl paraben, propyl paraben and purified water.

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This device is unitary packaged and sterilized with irradiation
Cannot be re-sterilized;
Expiring delay: 5 years (maintain performance and sterile state);

14. SYMBOLS

	Manufacturer		European Representative
	Manufacturing date		Do not re-sterilize
	Consult the instructions for use		Do not use if packaging is damaged
	Sterilized by irradiation		Reference
	Storage between 10 to 35°C		Batch number
	Expiry date		

15. Administrative information

Date of modification of the leaflet: Avril 2021
Date of first affixing of the CE marking: 2018

If you have any further questions or need information, please contact us at the following address:

Jiangsu Province JianErKang
Medical Dressing Co., Ltd.
NO. 1, Jianerkang Road
Zhixi Industry Concentration Area, Zhixi Town
213251 Jintan Changzhou City, China

Imported by:

ASEPT INMED
9, Avenue Mercure 31130 Quint Fonsegrives France
Tél : +33(0)5 62 57 69 00 Fax : +33(0)5 62 57 69 01
www.aseptinmed.fr

Shanghai International Holding Corp. GmbH
Eiffelstrabe 80, 20537 Hamburg, Germany