Product Description

This stethoscope is intended for single-patient use in isolation room settings. The stethoscope will eliminate the potential risk of patient-to-patient cross-contamination introduced by reusable stethoscopes. The stethoscope is available in adult and pediatric sizes. This IFU applies to model numbers: SPS-YA1010 (adult), SPS-YA1100 (adult), SPS-YP1010 (pediatric) and SPS-YP1100 (pediatric).

Intended Use

3M Single-Patient Stethoscope is intended for use by healthcare professionals for medical diagnostic purposes only. It can be used for auscultation of heart, lung, and other body sounds.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Precautions

- To avoid cross contamination between patients, do not use on more than one patient.
- To prevent ear canal injury, be sure the eartips are firmly attached to the headset before inserting them into your ears.
- Product contains small parts. To prevent choking hazard for children, keep out of reach of children.


Notice

To prevent stethoscope damage:

- Do not immerse the stethoscope in any liquid
- Do not subject the stethoscope to any sterilization process
- Avoid storing the stethoscope in extreme heat

Instructions for Use

Insert eartips in a forward direction.



The stethoscope headset is angled to complement the anatomy of the typical ear canal and is designed to provide a comfortable, acoustically sealed fit. The eartips should point in a forward direction as you insert them into your ears.

Listen to Low and High Frequency Sounds with the Pressure Sensitive Diaphragm

The chestpiece is equipped with a pressure sensitive diaphragm that enables you to emphasize either low or high frequency sounds by simply adjusting the pressure applied to the patient.

Low frequencies: apply light pressure to the chestpiece

High frequencies: apply firm pressure to the chestpiece

Cleaning

Eartips

- The stethoscope eartips should be cleaned between users with a 70% isopropyl alcohol wipe or a disposable wipe with soap and water. After cleaning, be sure the eartips are firmly attached to the headset before inserting them into your ears.

Stethoscope

- The stethoscope can be cleaned with 70% isopropyl alcohol wipe, a disposable wipe with soap and water, 2% bleach solution, or a quaternary disinfectant wipe.

Disposal: Dispose of contents/container in accordance with the local/regional/national/international regulations.

Explanation of Symbols

Symbol Title	Symbol	Description and Reference
Consult instructions for use		Indicates the need for the user to consult the instructions for use. Source: ISO 15223, 5.4.3
Natural rubber latex is not present		Indicates natural rubber or dry natural rubber latex is not present as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Authorized Representative in European Community		Indicates the authorized representative in the European Community. ISO 15223, 5.1.2
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. ISO 15223, 5.1.6
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. ISO 15223, 5.1.5
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. Source: ISO 15223, 5.1.1
Medical Device		Indicates the item is a medical device.
Unique device identifier		Indicates bar code to scan product information into patient electronic health record.
Importer		Indicates the entity importing the medical device into the EU.
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
CE Mark		Indicates conformity to European Union Medical Device Regulation or Directive.
MR Safe		An item that poses no known hazards resulting from exposure to any MR environment. MR Safe items are composed of materials that are electrically nonconductive, nonmetallic, and nonmagnetic. Source: ASTM F2503-13 Fig. 4
Intended for pediatric use		Indicates device is for use on pediatric patients

Made in the USA with Globally Sourced Materials

3M Company
2510 Conway Ave., St. Paul, MN 55144 USA
1-800-228-3957 (USA only)
3M.com

EC REP 3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str. 1
41453 Neuss, Germany

3M is a trademark of 3M.
Used under license in Canada.
©2019, 3M. All rights reserved.

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