3M™ Red Dot™
Monitoring Electrode with Foam Tape and Sticky Gel 2570 Series

Intended Use and Product Description:

The information provided in this insert applies to the 3M™ Red Dot™ Monitoring Electrode with Foam Tape and Sticky Gel 2570 Series (2570-3, and 2570-5) with snap connector and sticky gel. The 3M™ Red Dot™ Monitoring Electrodes are disposable with a silver/silver chloride sensing element and are designed for ECG monitoring. These electrodes have been tested for up to 5 days wear and are radiolucent.

Caution:

- The 3M™ Red Dot™ Monitoring Electrode with Foam Tape and Sticky Gel 2570 Series were determined to be MR-Conditional via non-clinical testing. For information on the use of 3M™ Red Dot™ monitoring electrode 2570 series during an MRI procedure please see below.
- U.S.A.: Federal law restricts this device to sale by or on the order of a physician.

Precautions:

- To prevent dryout, electrodes should be used within 45 days of opening the pouch. The unused electrodes should remain stored in the opened pouch.
- Do not use an electrode if the gel is dry or if the gel column is dislodged from the electrode.
- For proper skin management and to minimize skin irritation:
 - o Avoid placing an electrode on an irritated skin site.
 - o Do not abrade a skin site more than one time.
 - o Avoid removing electrodes frequently and/or reapplying to the same skin site.
 - o Avoid placing electrodes on skin wet from an alcohol wipe (dry thoroughly).
- Assess electrode sites periodically.
- During surgical procedures, place ECG electrodes as far as possible from the electrosurgical site to minimize RF current flow through the electrode site. Otherwise burns could result at the ECG electrode site.

Instructions For Use:

- Select and prepare an electrode site according to your health care facility's protocol for ECG monitoring or diagnostic procedures.
- The electrode site should be clean, dry and free of skin oil before application of the electrode to allow optimum electrode adhesion to skin.
- Application sites may be clipped.
- For Tape Preparation, the electrode should be part of the skin preparation.
- For Tape Preparation, the electrode should be part of the skin preparation.
- Attach the electrode to the electrode.
- Place the center of the electrode over the prepared skin site.
- Smooth the electrode down using a gloved motion. Avoid pressing on the center of the electrode.
- To minimize skin trauma during electrode removal, pull up slowly and gently from the edge of the electrode.

2570 MRI Information - MR Conditional

The 3M™ Red Dot™ Monitoring Electrode with Foam Tape and Sticky Gel 2570 Series was determined to be MR-Conditional via non-clinical testing.
A patient with this device can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MR-Related Heating

In non-clinical testing, the 3M™ Red Dot™ monitoring electrode 2570 series produced the following temperature rise during MRI performed for 15-min of scanning (i.e., per pulse sequence) in the 3-Tesla 0-Tesla/128-MHz, Excite, HDX, Software 14X.M3, General Electric Healthcare, Milwaukee, WI) MR system: Highest temperature change +1.5°C

Therefore, the MRI-related heating experiments for the 3M™ Red Dot™ monitoring electrode 2570 series, at 3-Tesla using a TransMx/veave RF body coil at an MR system reported whole body averaged SAR of 2.9-W/kg (i.e., associated with a calorimetry measured whole body averaged value of 2.7-W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.5°C.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the 3M™ Red Dot™ monitoring electrode 2570 series. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 2-mm relative to the size and shape of the 3M™ Red Dot™ monitoring electrode 2570 series.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
2570 Signal Body Size	1,483-mm²	130-mm²	1,822-mm²	273-mm²
Patient Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Important Note: This ECG Electrode that underwent MRI testing was not connected to an electrocardiograph (ECG) lead and/or an MR conditional monitoring system. Therefore, the safety of using this ECG Electrode attached to an ECG lead and/or MR conditional monitoring system is unknown. Furthermore, an MR conditional monitoring system typically requires use only with ECG electrodes specifically designed and tested with that specific equipment. It is unknown if this ECG Electrode is acceptable for use with an MR conditional ECG monitoring system designed to be used to monitor a patient undergoing an MRI procedure.

Shelf Life:


For shelf life, refer to the expiration date that is printed on each package.
Thank you for choosing 3M™ Red Dot™ monitoring electrodes. If you are not completely satisfied or if you have questions or comments, in the U.S.A., please call the 3M™ Health Care Customer Helpline at 1-800-228-3957, in Canada 1-800-364-3577. Outside of the U.S.A. and Canada, please contact your 3M™ Distributor or your 3M™ Sales Representative.

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3M and Red Dot are trademarks of 3M
© 2016, 3M. All rights reserved.