3M™ Red Dot™ ECG Monitoring Electrode with Foam Tape and Sticky Gel 2560 and 9643

Intended Use and Product Description:

The information provided in this insert applies to the 3M™ Red Dot™ ECG Monitoring Electrode with Foam Tape and Sticky Gel 2560 and 9643 with snap connector and sticky gel. The 3M™ Red Dot™ Monitoring Electrodes are disposable with a silver/silver chloride sensing element and are designed for ECG monitoring. These electrodes have been tested for up to 5 days wear.

Warning:

The 3M™ Red Dot™ monitoring electrodes 2560 and 9643 with a metal snap connector have not been tested for use during Magnetic Resonance Imaging (MRI) procedures. The use of conductive patient-connected devices and patient lead wires/electrodes in MRI procedures may result in serious patient burns. If a physician orders either of these devices to be used during MRI procedures, the physician should not allow loops in patient lead wires or allow the lead wires to come in contact with the patient's skin.

Caution:

U.S.A. federal law restricts this device to sale by or on the order of the physician.

Precautions:

- To prevent dryout, electrodes should be used within 45 days of opening the package. The unused electrodes should remain stored in the opened pouch.
- Do not use an electrode if the gel is dry or if the gel column is dislodged from the electrode.
- For proper skin management and to minimize skin irritation:
 - Avoid placing an electrode on an irritated skin site.
 - Do not abrade a skin site more than one time.
 - Avoid removing electrodes frequently and/or reapplying to the same skin site.
 - Avoid placing electrodes on skin still wet from an alcohol wipe (dry thoroughly).
 - Assess electrode sites periodically.
- During surgical procedures, place ECG electrodes as far as possible from the electrocauteral site to minimize RF current flow through the electrode site. Otherwise burns could result at the ECG electrode site.

Instructions For Use:

- Select and prepare an electrode site according to your health care facility's protocol for ECG monitoring or diagnostic procedures.
- The electrode site should be clean, dry and free of skin oil before application of the electrode to allow optimal electrode adhesion to skin.
- Application sites with heavy or curly hair that cannot be parted should be clipped.
- For best trace quality, skin abrasion should be part of the skin preparation. Use 3M™ Red Dot™ Trace Prep as part of the skin preparation.
- Attach the lead wire to the electrode.
- Remove electrode from the liner. Avoid touching the adhesive surface.
- Place the center of the electrode over the prepared skin site.
- Smooth the electrode down using a circular motion. Avoid pressing on the center of the electrode.
- To minimize skin trauma during electrode removal, pull up slowly and gently from the edge of the electrode.

Shelf Life:

For shelf life, refer to the expiration date that is printed on each package. Thank you for choosing 3M™ Red Dot™ monitoring electrodes. If you are not completely satisfied or if you have questions or comments, in the U.S.A., please call the 3M™ Health Care Customer Helpline at 1-800-228-3957, in Canada 1-800-364-3577. Outside of the U.S.A. and Canada, please contact your 3M™ Distributor or your 3M™ Sales Representative.

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