Periform[®]+

Instruction Manual – EN

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Manual de Instrucciones – ES

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Manuale d'Istruzioni – IT


Handleiding – NL

Bruksanvisning – SV

Brukerveiledning – NO

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inner strength, outer confidence

Specifications

Materials:

Probe body: High Impact Polystyrene

Electrodes: Medical Grade Stainless Steel

Cables: PVC over copper tinsel conductors

Weight: 19g including cables

Cable Length: 340mm

Dimensions: Length 94mm Width: 34.7mm

Connectors: 2 x 2mm diameter sockets

Electrode Surface Area: 4.9cm² per electrode

'In Use' product lifetime: 6 months

Contents

2 - 3	Instructions for Use - EN
4 - 5	Conseils d'utilisation - FR
6 - 7	Gebrauchsanweisung - DE
8 - 9	Instrucciones de Uso - ES
10 - 11	Instruções de utilização - PT
12 - 13	Istruzioni d'uso - IT
14 - 15	Gebruiksaanwijzing - NL
16 - 17	Bruksanvisning - SV
18 - 19	Brukerveiledning - NO
20 - 21	Käyttöohjeet - FI
22 - 23	使用说明 - ZH

Indications for Use:

This device should only be used with Neen Pericalm unit in the treatment of 'Stress', 'Urge' and 'Mixed' incontinence.

FOR SINGLE PATIENT USE – May be re-used by the same patient.

Use as directed by your medical adviser

Contraindications & Precautions:

DO NOT USE:

With electrical stimulation if fitted with a cardiac pacemaker • If pregnant
• If the package seal is not intact • If an infection or tissue damage or recent surgical scar is present • During menstruation • With electrical stimulation in patients with a history of carcinoma at the site of stimulation • If an irritation or discomfort occurs and seek medical guidance • If nursing. Do not use with children.

CAUTION IS ADVISED:

When using electrical stimulation on patients with altered or reduced tissue sensation • During electrical stimulation if the user has a history of epilepsy - seek informed medical guidance. If in doubt contact the equipment manufacturer or your medical adviser.

Cleaning:

- Wash PERIFORM⁺ in warm soapy water.
- Rinse in clean water and dry.
- Repeat after each use and dry thoroughly before storing.
- DO NOT use boiling water.
- Avoid immersion of the lead wires in water when cleaning.
- Do not autoclave.

Disposal:


Clean thoroughly and dispose of according to local regulations.

Contains nickel:

May cause irritation or allergic reaction following prolonged contact.

Always follow the user instructions.

Figure 1

Instructions for Use:

IMPORTANT:

- Before first and subsequent use wash and dry as directed.
- If using with Contraction Indicator connect together as shown in Figure 1.
- Insert into vagina in 'East/West' position, i.e. electrode surface facing towards hips; the external 'flange' should sit comfortably between the labia.
- Ensure equipment is turned off before connecting to probe cables.
- Use as directed by your medical adviser. After use, turn off equipment, disconnect cables. Do not pull on cables to remove probe.
- Wash, dry and store as instructed.


Contraction Indicator:

When fitted to the PERIFORM[®]+ the Contraction Indicator provides valuable evidence of the performance of the pelvic floor muscles. It may be used on its own or when using stimulation and biofeedback equipment. When pelvic floor muscles contract correctly the indicator will move downwards (✓).

Moving upwards indicates incorrect contraction (X) as shown in Figure 2.

Functions and benefits of the Indicator are:

- Quantifies the strength and duration of voluntary and involuntary (stimulated) contractions. N.B. Voluntary contractions can be used to reinforce stimulated contractions.
- Indicates need to use stimulation to strengthen a weak pelvic floor – i.e. little or no indicator movement.
- Teaches bracing techniques to help avoid leakage when coughing / sneezing / straining etc.
- Demonstrates the success of a pelvic floor exercise program.

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
Notes

Notes

For more information on this and other products in the NEEN® range please visit:
www.neenpelvichealth.com

inner strength, outer confidence

Serious events experienced from using this product should be reported immediately to Performance Health and to the local Competent Authority.

These user instructions are available in large print, please contact Performance Health Customer Services for a copy. Please retain these instructions for future reference.

Performance Health International Ltd.
 Nunn Brook Road,
 Huthwaite, Sutton-in-Ashfield,
 Nottinghamshire, NG17 2HU, UK.

UK:
 Tel: 03448 730 035
 Fax: 03448 730 100
www.performancehealth.co.uk
 International:
 Tel: +44 1623 448 706
 Fax: +44 1623 448 784

Distributed in Australia by:
 Performance Health ANZ
 Unit 3, 3 Basalt Road,
 Pemulvuy NSW, Australia 2145.
 Tel: 1300 473 422
 Fax: 1300 766 473
www.performhealth.com.au

Distributed in France by:
 Performance Health France
 13 rue André Pingat,
 CS 10045
 51724 Reims Cedex
 France
 Tel: 00(33) 03 10 00 79 30
 Fax: 00(33) 03 10 00 79 49
www.performancehealth.fr

CEpartner4U
 Esdoornlaan 13
 3951 DB Maarn
 The Netherlands
www.cepartner4u.com

71524 1s7

NEEN® & PERIFORM® are registered trademarks of Performance Health International Limited. 2021