
## MED-GEL

|    |                              |
|----|------------------------------|
| en | INSTRUCTIONS FOR USE         |
| de | GEBRAUCHSANWEISUNG           |
| fr | NOTICE D'UTILISATION         |
| nl | GEbruikSAANWIJZING           |
| it | ISTRUZIONI PER L'USO         |
| es | INSTRUCCIONES DE USO         |
| pt | INSTRUÇÕES PARA USO          |
| sv | INSTRUKTIONER FÖR ANVÄNDNING |
| fi | KÄYTTÖOHJEET                 |
| da | INSTRUKTION                  |
| no | BRUKSANVISNING               |

15 g

2 Kg

5 °C

**Medline International France SAS,**

5 rue Charles Lindbergh, 44110 Châteaubriant, France.

Distributed in Australia by: Medline International Two Australia Pty. Ltd., Sydney, Australia

[www.medline.com](http://www.medline.com)

RB17HBN 2017-02

### Description:

MED-GEL is a granular substance (Sodium polyacrylate, cross linked) which jellifies aspirated fluids in the Operating Theatre and Intensive Care Unit to facilitate safe disposal of waste.

### Composition/Ingredient:

Sodium Polyacrylate (95%), Water (5%).

### Recommended Dosage:

RE2015: MED-GEL is packed in an individual 15gr water-soluble bag. 1 bag = 1 dose = 15gr MED-GEL

RE2050: MED-GEL is packed in a 2kg container which holds 2 x1kg bag and a measuring cup. Use the measuring cup to pour MED-GEL to jellify the fluids. 1 full measuring cup = 1 dose = 15gr MED-GEL

| Volume of fluids to be jellified | Quantity of MED-GEL required | Equivalence RE2015 | Equivalence RE2050 |
|----------------------------------|------------------------------|--------------------|--------------------|
| 1.0 L                            | 15 gr of MED-GEL             | 1 dose of 15 gr    | 1 cup of 15 gr     |
| 1.5 L                            | 30 gr of MED-GEL             | 2 doses of 15 gr   | 2 cups of 15 gr    |
| 3.0 L                            | 45 gr of MED-GEL             | 3 doses of 15 gr   | 3 cups of 15 gr    |

The fluids are jellified after approximately 10 - 15 minutes.

Dosage is dependent on the composition of the fluids. A high percentage of electrolytes in the solution typically requires a higher dosage. Gallbladder fluids can partially reduce the effectiveness (fat).

### Recommended Application: with disposable suction liner and canister system

MED-GEL is added in the appropriate quantity to the fluids in the disposable suction liner or canister. It is recommended to add MED-GEL (via the accessory port) into the suction liner after suctioning and prior switching off the vacuum source.

- In single configuration: MED-GEL may be added into the liner prior to suctioning. If premature shutoff occurs it will not be reversible and will require the immediate change of the liner.
- In tandem configuration: Do not put MED-GEL into the liner prior to suctioning since this may reduce suction effectiveness. Add MED-GEL after suctioning and prior switching off the vacuum source.

### First aid measures:

- Skin Exposure: In case of contact, immediately wash skin with soap and copious amounts of water. If irritation persists, call a physician.
- Eye exposure: In case of contact with eyes, flush with copious amounts of water for at least 15 minutes. Assure adequate flushing by separating the eyelids with fingers. If irritation persists, call a physician.
- Inhalation Exposure: If inhaled, away from the area of exposure. If necessary, seek medical attention.
- Oral Exposure: If swallowed, irrigate the mouth with water and seek medical attention.

### Exposure controls/Personal protection:

- Respiratory protection: Use respiratory protection when dust is formed. Ensure sufficient workplace ventilation.
- Eye protection: Protective goggles
- Hand protection: Not applicable

**Fire-fighting measures:**

- Extinguishing media: Suitable: water spray, dry chemical, carbon dioxide or appropriate foam.
- Firefighting protective equipment: wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.

**Toxicity, Dangers:**

- Acute oral toxicity: LD50, dose > 5000 mg/kg, limit test.
- Skin irritation: No irritation.
- Eye irritation: Very slight eye irritation. Particle effect.
- Sensitization: No sensitization.
- Acute dermal toxicity: LD50, dose > 2000 mg/kg, limit test.

**Disposable Considerations:**

- To be stored under dry conditions.
- Do not dispose of jellified fluids via the sewage system. Comply with facility's protocols and local, state and federal requirements for disposal (COSHH regulations for UK only).
- In regards to the present knowledge the product is a non-dangerous substance in the definition of the European Regulation N°1272/2008 – Classification, Labeling and Packaging of substances and mixtures.

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