Flexible Enema Administration Set, with soap, underpad, and prelubricated tip.

Instruction for use

- de** **Flexibles Einlauf-Verabreichungs-Set**
Gebrauchsanweisung
- fr** **Ensemble pour lavement**
Notice d'utilisation
- nl** **Flexibele klysματοedieningsset**
Gebruiksaanwijzing
- it** **Set di somministrazione clistere flessibile**
Istruzioni per l'uso
- es** **Set flexible para administración de enemas**
Instrucciones de uso
- pt** **Conjunto flexível de administração de enema**
Instruções para Uso
- sv** **Flexibelt administreringsset för lavemang**
Instruktioner för användning
- da** **Flexibelt lavementsæt**
Instruktion
- no** **Flexibelt klysteradministrasjonssett**
Bruksanvisning
- pl** **Elastyczny zestaw do wykonywania lewatywy**
Instrukcja obsługi
- ro** **Set flexibil de administrare a climei**
Instrucțiuni de folosire

Conod Medical Co., Ltd,
No. 11 Hongfeng Road, Baimao Industrial Park, Guli Town,
215532 Changshu City, Jiangsu Province, P.R. China

Shanghai International Holding Corp. GmbH,
Eiffestraße 80, 20537 Hamburg, Germany

Medline International B.V.,
Nieuwe Stationsstraat 10, Arnhem, 6811 KS, Netherlands

CH Rep and importer:
The Tao of Excellence GmbH, Vorstadt 26,
CH-8200 Schaffhausen, Switzerland

RA23CSI 2023-01

REF DYNDE70100A

Instruction
for use

en Flexible Enema Administration Set, with soap, underpad, and prelubricated tip

Contents:

1500 ml Enema Bag w/preattached Tube Lubricated Tip
Tip Protector
Tubing Clamp
Underpad
Soap Packet

Directions For Use:

Slide tubing clamp to «closed» position.
Open top of bag and add desired quantity of solution.
Close top of bag with thumb and forefinger to cause the sides of the valve closure to stick together.
Suspend bag in preparation for administration.
Position moisture proof drape under patient.
Remove tip protector from tubing prior to use.
Proceed with administration using approved technique.
Drain unused solution from bag and discard.
Wash hands before and after procedure.

Caution:

Water or enema solution at high temperature is forbidden to avoid damage to the anus and intestinal mucosa.
Pre-checking the rectal tube of bag/bucket before use whether open tip or side hole of tube is blocked by lubricant.
Remove tip protector from tubing prior to use.
This device should only be used by trained healthcare professionals.

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en If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established. Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste. Reuse of single use devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function; which may lead to injury, illness and/ or death of a patient.

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