Lomatuell® Pro

Gelbildendes Kontaktnetz
Interface gélifiante
Gel-forming contact net
Malla de contacto gelificante
Rede de contacto formadora de gel
Rete di contatto gelificante
Gelvormend contactnetwerk
Gelbildende kontaktnät
Geelyityvä haavaverkko
Kontaktní síťka vytvářející gel
Kontaktni sićka tvoriaca gél
Opatrunek sićkowy z żelującą warstwą kontaktową
Kontaktna mreža, ki tvori gel
Kontaktna gelirajuća mrežica
Πλέγμα επαφής για σχηματισμό γέλης
Plasa de contact cu formare de gel
Geeli moodustav kontakt-võrkside
Gelį formuojantis kontaktinis tinklielis

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en Instructions for Use

Product description and features

Lomatuell® Pro is a non-adhesive, flexible wound dressing with an open net structure and good exudate draining into the secondary dressing (wound contact layer). The product consists of a synthetic tulle saturated with a polymer matrix that contains lipids and hydrocolloids. On contact with fluids, the hydrocolloid particles form a gel, which forms an amphiphilic (hydrophilic-lipophilic) layer with the petroleum jelly. Adherence to the wound is largely prevented, wound healing is supported and traumatic dressing changes are facilitated.

Lomatuell® Pro can be used in combination with negative pressure therapy as wound contact layer* (superficial wounds) and it can be cut to size for different wound sizes and areas of the body.

*based on in-vitro studies

Product composition

Petroleum jelly, medical-grade white oil, polyester (polyethylene terephthalate), sodium carboxymethyl cellulose, synthetic rubber, antioxidant pentaerythritol tetrakis(3-(3,5-di-tert-butyl-4-hydroxyphenyl)propionate)

Indications

Lomatuell® Pro is suitable for the following indications :

Chronic wounds :

- lower leg ulcers,
- diabetic foot ulcers,
- pressure ulcers.

Acute wounds:

- lacerations, cuts and abrasions,
- second-degree burns.

Post-surgical wounds:

- surgical wound healing by secondary intention,
- skin graft donor sites.

Contraindications

Known allergy and/or hypersensitivity to any of the product components

Side effects

In very rare cases skin irritations and/or allergies may occur.

Warnings and precautions

Lomatuell® Pro is sterile as long as the packaging is unopened and undamaged. Do not use if package is damaged or accidentally opened. The product is intended for single use and must not be resterilised.

Instructions for use

The product is applied by healthcare professionals.

1. Preparation :

Before applying Lomatuell® Pro, carefully cleanse the wound by irrigating it (e.g. with physiological saline or Ringer's solution).

2. Application of the wound dressing:

Before opening, check the sterile barrier for obvious damage

- Open the sterile packaging at the designated place and remove the product in a sterile manner.

- Remove the cover film, which is applied to both sides

- Apply Lomatuell® Pro so it extends past the wound by around 1 cm. Where needed, the product can be cut to size with sterile scissors to adapt the dimensions of the wound dressing to the wound.

- For highly exuding wounds, use a highly absorbent wound dressing as secondary dressing (e.g. Vliwazel® absorbent dressing), for moderately to slightly exuding wounds, a gauze or nonwoven dressing pad (e.g. Gazin® gauze swabs or Vliwasoft® nonwoven wound dressing) is sufficient. To promote a moist, warm wound environment and protect the wound from external germs and bacteria, the wound dressing may be retained with a PU film (e.g. Suprasorb® F film dressing). In case of increased risk of wound infection or infected wounds, dressing retention should be performed with air-permeable materials (e.g. Curafix® H dressing retention sheet). In the case of infected wounds, the established procedures for treatment of infections must be carried out in parallel.

3. Dressing change:

The recommended dressing change interval is normally 2–4 days. Depending on the type of wound and its healing progress, the dressing can be left on the wound for a shorter or longer period (up to a maximum of 7 days) – in consultation with the responsible medical professional.

4. Storage and transport conditions:

Store Lomatuell® Pro horizontally.

5. Disposal:

In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance – AVV). Packaging labelled as recyclable should be taken into the respective national recycling systems.

General instructions

In case of a serious incident, contact the manufacturer and the responsible health authorities.

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