Rosidal® sys simple

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en Instructions for Use

Product description and performance characteristics

Compression system for venous leg ulcer therapy
Rosidal® sys simple is available in two versions with components of differing quantities and sizes to individually adapt the therapy to the patient.

Product composition

147647 Rosidal® sys simple 8/10

1 x Rosidal® soft (foam bandage, white),
10 cm x 2.5 m x 0.3 cm
100% polyurethane

1 x Rosidal® K (short stretch bandage, cream-coloured),
8 cm x 5 m
100% cotton

1 x Rosidal® K (short stretch bandage, cream-coloured),
10 cm x 5 m
100% cotton

1 x Mollelast® haft (cohesive conforming bandage, white),
8 cm x 4 m
Backing: 66% viscose, 34% polyamide,
natural rubber adhesive

147648 Rosidal® sys simple 10/10

2 x Rosidal® soft (foam bandage, white),
10 cm x 2 m x 0.2 cm
100% polyurethane

2 x Rosidal® K (short stretch bandage, cream-coloured),
10 cm x 5 m
100% cotton

1 x Mollelast® haft (cohesive conforming bandage, white),
8 cm x 4 m
Backing: 66% viscose, 34% polyamide, natural rubber
adhesive

Intended purpose

For venous leg ulcer compression therapy

Indications

Venous leg ulcers and their comorbidities, as well as chronic oedema (due to venous insufficiency)

Contraindications

Advanced peripheral arterial occlusive disease, decompensated heart failure, septic phlebitis, phlegmasia coerulea dolens, sensation disorders of the skin, known allergies and/or hypersensitivity to product components.

Side effects

In very rare cases skin irritations and / or allergies may occur.

Warnings and precautions

- Incorrect bandaging can result in skin necrosis and pressure damage to nerves.
- Slight bluish tinge to the toes is normal, this must disappear as soon as the patient starts to walk.
- If the symptoms (pain and bluish tinge) do not disappear, the bandage must be removed and applied again with less compression.
- The patient must inform their doctor if sharply increasing pain occurs while wearing the bandage.
- Use with an Ankle Brachial Pressure Index (ABPI) of over 0.8, or with an ABPI of 0.5-0.8 – under medical supervision only.
- In the case of diabetes with advanced microangiopathy – use only under medical supervision.

Application instructions

General:

- For reuse on a single patient.
- The compression bandage can stay on the leg for up to 7 days, depending on the wound status.
- The application technique is chosen in consultation with the treating physician.
- Always apply Rosidal® K with full tension (max. stretch).
- It is recommended to apply the bandage in the morning before the patient gets up or after the patient has rested the leg in a raised position for a prolonged period and to then encourage the patient to walk.
- The width of the bandage should correspond to the diameter of the extremity to be bandaged.
- The bandage must be applied with tension around the leg in such a way that the pressure from the bandage decreases from the foot to the knee and, where applicable, to the gluteal fold.
- For more skin protection, a tubular bandage (e.g. tg® Tubular Bandage) can be applied additionally under the compression bandage.
- Tubular bandage and adhesive strips are not part of this set.

Preparation:

- The usual diagnostic examinations to determine the vascular status must be carried out before applying the bandage.
- Before applying the bandage, clean/debride any wound that is present and cover with a wound dressing appropriate for the wound status.

Application instructions:

- Training by healthcare professionals on the correct application is required for the untrained user.
- The bandage is unrolled with the rolled-up portion (roll) on top. Keep the roll close to the leg.
- To achieve maximum joint movement, it is recommended to apply the bandage with the foot at a right angle (90°) to the leg.
- Use adhesive strips to fix the bandage in place. Do not fix in place with bandage clips because these can cause injury.

Bandaging instructions (modified after Sigg) (see Figures)

- Start applying Rosidal® soft (padding layer) on the top of the foot going around the inside of the instep including the base of the toes.
- The wraps of the bandage should always overlap by approx. ½ - 2/3.
- Spiral wrap to the heel including the heel.
- Return down around the ankle to secure the heel (towards the sole).
- To finish securing the heel, wrap above the ankle (towards the Achilles tendon).
- Continue to spiral wrap to the base of the calf muscle.
- Depending on the version of the set, wrap the calf using a figure of eight pattern with the second Rosidal® soft.
- Secure the end of the bandage below the kneecap by wrapping the Mollelast® haft once around.
- Start applying the compression bandage Rosidal® K (8 cm or 10 cm wide, depending on the version of the set) at the top of the foot going around the inside of the instep including the base of the toes.

- The wraps of the bandage should always overlap by approx. ½ - 2/3. Always apply Rosidal® K with full tension (max. stretch).
- Spiral wrap to the heel including the heel.
- Return down around the ankle to secure the heel (towards the sole).
- To finish securing the heel, wrap above the ankle (towards the Achilles tendon) and continue to spiral wrap to the base of the calf muscle.
- Wrap the calf using a figure of eight pattern with Rosidal® K (10 cm wide).
- Secure the end of the bandage below the kneecap with an adhesive strip (e.g. Porofix®).
- Secure the tubular bandage below the knee and on the forefoot by wrapping Mollelast® haft once around.
- The finished bandage.

Care recommendations

Rosidal® K can be washed up to 50 times at 95°C (standard commercially available mild detergent/delicates cycle).
Dryer-safe at max. 55°C, set iron to low heat.
Rosidal® soft can be washed up to 50 times at 95°C (standard commercially available mild detergent/delicates cycle). Not dryer-safe, do not iron.

The maximum stretch of the bandage may increase with washing. Please bear this in mind when applying the bandage.

Disposal

In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Recyclable packaging should be recycled according to the relevant national recycling systems.

General instructions

In case of a serious incident, contact the manufacturer and the responsible health authorities.

Contains or presence of natural rubber latex.