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Ferricyanure de potassium 2,2 % p/p

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
Références

1. American Diabetes Association, Standards of Medical Care in Diabetes, Diabetes Care, 2016:39 (Supplement 1) S39-S46.
2. Atkins SH, Dasmahapatra A, Jaker MA, Chorost MI, Reddy S: Fingertstick glucose determination in shock. Ann Int Med 114: 1020-1024. 1991
3. Sandler M, Low-Beer T: Misleading capillary glucose measurements. Practical Diabetes 7: 210. 1990
4. Wickham NWR, Achar KN, Cove DH: Unrealiability of capillary blood glucose in peripheral vascular disease. Practical Diabetes 3: 100. 1986.

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Bandelette de test

Spengler SAS
30 rue Jean de Guirmand
13290 Aix-en-Provence - FRANCE

Tyson Bioresearch, Inc.
5F., No. 16,18,20,22, Kedong 3rd Rd.,
ZhuNan Township, Miaoli County 35053,
Taiwan

Medical Device Safety Service GmbH (MDSS)
Schiffgraben 41, D-30175 Hannover, Germany

Centre de support technique client :
+33(0)4 42 90 31 31 (choix 83)

Spengler JT100 Test Strip Package Insert

For Testing Glucose in Capillary and Venous Whole Blood

IMPORTANT: The JT100 Blood Glucose Monitoring System is comprised of the JT100 Blood Glucose Meter and JT100 Blood Glucose Test Strips. Please read this information and your JT100 user's manual before using JT100 Test Strips. Call Spengler Customer Service immediately if these instructional materials or your meter results seem unclear.

Intended Use

The JT100 Blood Glucose Test Strips are used with the JT100 Blood Glucose Meter for quantitatively measuring glucose (sugar) in whole blood obtained from the fingertip, palm and forearm, or venous whole blood. The JT100 Blood Glucose Monitoring System is intended for use in the home and in professional settings to monitor blood glucose levels. **The JT100 Test Strips are for testing outside the body (in vitro diagnostic use).**

Test Principle

On each test strip there is a test area containing reaction chemicals. When blood is applied to this area, a chemical reaction takes place, then a transient electrical current is formed. The blood glucose concentration is calculated based on the electrical current detected by the meter, then the result is shown on the display. The test measures glucose from 10mg/dL (0.56mmol/L) to 630 mg/dL(35mmol/L). The JT100 Test Strip is calibrated to display the equivalent of plasma glucose values to allow easy comparison of results with laboratory methods.

Precautions

- **For in vitro diagnostic use (outside the body) only.**
- Do not reuse test strips.
- Do not use test strips beyond the expiration date printed on the package since this may cause inaccurate results.
- Discard the vial and any unused test strips 12 months after you first open it. Constant exposure test strips to air may destroy chemicals in the test strip. This damage can cause incorrect readings.
- The sample must only be applied to the tip of the test strip. Do not apply blood or control solution to the top, left or right sides of the test strip. This may cause inaccurate test results.
- Use only the JT100 Blood Glucose Meter with the JT100 Test Strips.
- Keep the test strip vial away from children. The cap is a choking hazard. The cap or vial contains drying agents to protect the test strips. Drying agents may be harmful if inhaled or swallowed and may cause skin or eye irritation.
- **Any change in medication based on the JT100 blood glucose test results without the consent and advice of a physician or healthcare professional is not recommended. Consult your physician or healthcare professional before making changes in your treatment plan based on your blood glucose test results.**
- Severe dehydration and excessive water loss may cause false low results. If you believe you are suffering from severe dehydration, consult your physician immediately.
- All devices contaminated with blood should be disposed of properly.
- Healthcare professionals should follow their institution's infection control protocols.
- If you get JT100 Control Solution test results that fall outside the range printed on the test strip vial, repeat the test with a new test strip. If the test result is still outside the range printed on the test strip vial, the JT100 may not be working properly. DO NOT use the system to test your blood until you get a control solution test result within the range printed on the test strip vial label.

Storage and Handling

- Store at room temperature between 2-32°C (36-90°F). Do not refrigerate or freeze. Avoid exposing test strips to extreme temperatures.
- Store away from direct sunlight and heat.
- Use each strip immediately after removing it from the vial.
- Store your test strips in their original vial only. The cap or vial contains drying agents to protect the test strips from damp. Do not transfer test strips to a new vial or any other container.
- **After removing a test strip from the vial, replace the vial cap immediately and close it tightly.**
- Do not bend, cut, or alter the Test Strip in any way.
- With clean, dry hands you may gently touch the test strip anywhere when removing it from the vial or inserting it into the meter.

How to do the Test

See your User's Manual for a step-by-step guide on how to do the test.

1. Select the Puncture Site to Obtain Blood

You may obtain a blood drop from the fingertip, palm, or forearm. If using the forearm, choose a fleshy area away from bone, visible veins and hair. If using the palm, choose a fleshy area on the palm below your thumb or pinky finger. Select a spot with no visible veins and away from deep lines, which may cause your blood sample to smear. Wash your hands and the puncture site with warm, soapy water. Dry them thoroughly. If you use an alcohol swab, make sure the skin is dry before lancing. To help you get a large enough drop of blood, you may gently massage or apply warmth to the site to increase blood flow. Press and hold the sampler against your fingertip, palm or forearm for a few seconds, then press the release button.

Warning

Alternate site testing (AST) should ONLY be used in the following intervals:

- In a pre-meal or fasting state, (more than 2 hours since the last meal)
- Two hours or more after taking insulin.
- Two hours or more after exercise.

DO NOT use AST if:

- You think your blood glucose is low.
- You are unaware of hypoglycemia.
- Your AST results do not match the way you feel.
- You are testing for hyperglycemia.
- Your routine glucose results are often fluctuating.
- If you are pregnant.

All parts of the glucose monitoring system should be considered potentially infectious and are capable of transmitting blood-borne pathogens between each user. To assure that you are not placing yourself or other persons at risk, always remember:

- The meter and lancing device are for single person use.
Do NOT share them with anyone including other family members! Do not use on multiple persons!
- All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection.
- Refer to your user's manual section "Cleaning and Disinfection procedure" for how to clean your meter and lancing device.

2. Blood Glucose Testing

Insert the test strip into the JT100 blood glucose meter. The meter will turn on and detect automatically, please go to next step.

3. Use the Lancing Device to obtain the right size blood drop. Refer to your Lancing Device User's Instruction for how to lance your finger, forearm and palm. Then obtain a blood sample about the size of a pinhead. Gently touch only one edge of the Test Strip to the blood sample. You may fill the blood to the tip of the test strip. Do not apply blood to test strip except the tip of test site. Do not press the edge of the strip against the test site. When the strip is full of blood, the meter will "beep" or you will see the test strip chamber area completely filled.

4. Read Results

Read the test results in an average of 5 seconds. The result will be stored in the meter memory. Results are displayed as mg/dL (milligrams per deciliter of blood) or mmol/L (millimoles per liter of blood) depending on the units you chose when you set up your meter. To change the units displayed, see your User's Manual section "Setting Up the Meter."

Wash hands thoroughly with soap and water after handling the meter, lancing device, or test strips.

What Does The Results Mean?

Normal Glucose Values

The glycemic recommendations for non pregnant adults with diabetes :

- Preprandial capillary plasma glucose : 80–130mg/dL (4.4–7.2mmol/L)
- Peak postprandial capillary plasma glucose: 180mg/dL (10.0mmol/L)

For people with diabetes, consult your physician or healthcare professional for the target glucose values that are right for you.

Low Glucose Values

The JT100 Meter displays results between 10 and 630 mg/dL (0.56 and 35 mmol/L). If your test result is lower than 10, "Low" (LO) will appear on the meter display. This indicates severe low blood sugar (hypoglycemia). You should immediately treat low blood sugar as recommended by your healthcare professional.

High Glucose Values

If your test result is above 630 mg/dL (35 mmol/L), "High" (HI) will appear on the meter display screen. This indicates severe high blood sugar (hyperglycemia). You should immediately treat high blood sugar as recommended by your healthcare professional.

Unexpected Results

Low or high blood sugar readings can indicate a potentially serious medical condition. If your blood sugar is unusually low or high, or if you do not feel the way your results indicate, repeat the test with a new test strip. **If your reading is not consistent with your symptoms or if your blood glucose result is less than 60 mg/dL (3.3 mmol/L) or higher than 240 mg/dL (13.3 mmol/L), you should contact your healthcare professional and follow his or her treatment advice.**

Limitations

The JT100 Blood Glucose Test Strips give accurate results when the following limitations are observed:

1. Use only the JT100 Test Strip with the JT100 Meter.
2. Use only capillary capillary and venous whole blood from finger, palm and forearm. Do not use plasma or serum or blood samples from other than capillary or venous whole blood.
3. To minimize the effect of glycolysis, blood glucose determination with venous blood must be performed within 30 minutes of sample collection.
4. Fresh venous blood samples containing heparin is acceptable.
5. Do not use for testing neonatal blood.
6. The system is tested to accurately read the measurement of glucose in whole blood within the range of 10 to 630mg/dL (0.56-35 mmol/L).
7. This system can be used up to an altitude of 10,745 feet (3,275 meters).
8. The test strips are for single use only. Do not reuse test strips.

Additional information for Healthcare Professionals

9. JT100 Test Strip does not interfere with the hematocrit at a range of 0–70% of blood glucose. High (above 70%) hematocrit can cause false results. Patients with higher hematocrit values (above 70%) will usually test lower. Please ask your healthcare professional if you do not know your Hematocrit level.
10. Acetaminophen, uric acid, ascorbic acid (vitamin C), and other reducing substances when occurring in normal blood or normal therapeutic concentrations do not significantly affect results. However, abnormally high concentrations in blood may cause inaccurately high results. Lipemic samples; Cholesterol up to 1,200 mg/dL or triglycerides up to 1,525 mg/dL do not significantly affect test results. However, glucose values in specimens beyond these levels, should be interpreted with caution.
11. Antiglycolysis and anticoagulants in blood samples may affect the test results.
12. Severe dehydration, diabetic ketoacidosis, hyperosmolar non-ketotic state, hypotension, shock, or peripheral vascular disease 2,3,4 may cause inaccurate results.
13. Critically ill patients should not be tested with home use blood glucose monitors as it may lead to inaccurate results.
14. To minimize the effect of glycolysis, blood glucose determination with venous blood is suggested to be performed within 30 minutes of sample collection.
15. Fresh venous blood samples containing heparin is acceptable.

Summary of Substances with Interference

Interfering substance	The max toleration concentration in the evaluation
Acetaminophen	4.25 mg/dL
Ascorbate (Ascorbic acid)	3.0 mg/dL
Bilirubin	25.0 mg/dL
Cholesterol	1200.0 mg/dL
Creatinine	10.0 mg/dL
Dopamine	2.0 mg/dL
EDTA	200.0 mg/dL
Galactose	500.0 mg/dL
Genitiscic acid	2.5 mg/dL
Glutathione	100.0 µmol/L
Haemoglobin-HUMAN	3.0 g/dL
Heparin	5.0 IU/mL
Ibuprofen	500.0 µg/mL
Icodextrin	750.0 mg/dL
L-DOPA(L-3,4-dihydroxyphenylalanine)	5.0 µg/mL
Maltose	2575 mg/dL
Methyldopa	3.0 mg/dL
Pralidoxime iodide (PAM)	5.0 mg/dL
Salicylate	500.0 µg/mL
Tetracycline	5.0 mg/dL
Tolbutamide	100 mg/dL
Tolazamide	6.0 mg/dL
Triglycerides	1525.0 mg/dL
Urate(Uric acid)	8.0 mg/dL
Xylose	5.0 mg/dL
Sodium Carbonate	37.5 mEq/L
Mannitol	0.09 mg/dL
Sorbitol	0.09 mg/dL
Xylitol	0.09 mg/dL
Isomalt	0.09 mg/dL
Lactitol	0.09 mg/dL
Maltitol	0.09 mg/dL
Hydrogenated starch hydrolysates (HSH)	0.09 mg/dL

Quality Control

The control solutions are used to check the performance of the JT100 meter, Test Strips, and your testing technique. The system is performing correctly if the control solution test result falls within the specific control solution range listed on your Test Strip vial. A control solution test should be performed when

- You get your meter for the first time before testing your blood test.
- You open a new vial of test strips.
- You want to check the meter and strips.
- You suspect that the meter or test strips are not working properly.
- Your blood glucose test results are not consistent with how you feel.
- You think your test result does not accurate.
- You dropped your meter.
- Whenever your healthcare professional recommends it.

When a control solution test is done, you should get results within the expected range printed on the test strip vial. If control solution test results fall outside this range, repeat the test. Results that fall outside the range may be caused by:

- error in performing the test,
- expired or contaminated control solution/tests strips,
- test strip deterioration or meter malfunction.

If the problem persists, the system may not be working properly. DO NOT continue to use the system to test your blood glucose level. If you continue to get results that fall outside the specified range, please call Customer Service.

Note: JT100 Control Solutions are required but not supplied with Test Strips for checking the system. For order information, please contact Customer Service. There are three JT100 Control Solutions at different desired range of aqueous glucose (1, 2 or 3). Two sets of control solution should be used, a high level and a low level control solution of your choice. For more information about the Control Solution, please read your Control Solution package insert.

Performance Characteristics
Measuring Range
10-630 mg/dL (0.56-35 mmol/L)

System Accuracy
The accuracy of the JT100 System was compared with YSI 2300 Stat Plus glucose analyzer with fingertip testing done by persons with diabetes. ISO 15197: 2013 requirement: 95% of individual glucose results falling within ±15mg/dL at glucose concentration for samples <100mg/dL and within ±15% at glucose concentrations≥100mg/ dL. The total accuracy study was performed with:

- 100 people
- 3 lots of test strips
- the total measured values is 600
- 594 of 600 values (99.0%) are within the system accuracy criteria

Totally analysis system accuracy of lot1, lot 2 and lot 3.

Accuracy results for glucose concentration < 100mg/dL(5.55 mmol/L)		
Within ± 5mg/dL (0.27 mmol/L)	Within ± 10mg/dL (0.56 mmol/L)	Within ± 15mg/dL (0.83 mmol/L)
60,3 % (123/204)	91,2 % (186/204)	100,0 % (204/204)

Accuracy results for glucose concentration ≥100mg/dL (5.55 mmol/L)		
Within ± 5%	Within ± 10%	Within ± 15%
60,6 % (240/396)	87.1%(345/396)	98,5 % (390/396)

Accuracy results for glucose concentrations between 27 and 522 mg/dL		
Within ±15mg/dL(0.83 mmol/L) or ±15 %		
99,0 % (594/600)		

Consensus Error Grid of Three Lots						
Area	A+B	A	B	C	D	E
Points	600	600	0	0	0	0
CEG (%)	100 %	100 %	0 %	0 %	0 %	0 %

User Performance Evaluation
The lay user study of JT100 blood glucose monitoring system was assessed with:

- 1 lot of test strips
- 106 lay users
- 2 meters 'readings
- compared with YSI2300 glucose analyser

A study evaluating glucose values from fingertip capillary blood samples obtained by 106 lay persons showed the following results:

Accuracy results for glucose concentration < 100mg/dL(5.55 mmol/L)			
	Within ± 5 mg/dL (0,28 mmol/L)	Within ± 10 mg/dL (0,56 mmol/L)	Within ± 15 mg/dL (0,83 mmol/L)
Fingertip	50,0 % (6/12)	91,7 % (11/12)	100 % (12/12)
Palm	33,3 % (4/12)	83,3 % (10/12)	100 % (12/12)
Forearm	66,7 % (8/12)	75,0% (9/12)	100 % (12/12)

Accuracy results for glucose concentration ≥ 100mg/dL (5.55 mmol/L)			
	Within ±5%	Within ±10%	Within ±15%
Fingertip	47,9% (45/94)	85,1% (80/94)	96,8% (91/94)
Palm	46,8% (44/94)	85,1% (80/94)	98,9% (93/94)
Forearm	42,6% (40/94)	80,9% (76/94)	97,9% (92/94)

Accuracy results for glucose concentrations between 66 and 374 mg/dL	
	Within ± 15 mg/dL (0,83 mmol/L) or ±15 %
Fingertip	97,2%(103/106)
Palm	99,1% (105/106)
Forearm	98,1% (104/106)

Repeatability and Precision
Between-run precision of JT100 test strips was measured with venous blood samples in the laboratory. The results for total three hundred test strips from 3 lots (n=300) concluded the following repeatability and precision estimates.

Mean Blood Glucose (mg/dL)	36.9	83.5	122.9	231.9	331.3
SD (mg/dL) or CV (%)	2.2	2.8	2.4%	2.4%	2.3%

Intermediate Precision
Ten replicate assays drawn from 3 strip lots were run on JT100 Blood Glucose Meter each day for a total of 10 days. Control solutions at three concentration levels were used in the testing. The results concluded the following intermediate precision estimates.

Mean of control solution (mg/dL)	46.9	122.6	345.0
SD (mg/dL) or CV (%)	1.5	1.7%	1.8%

Chemical Composition
Glucose dehydrogenase (FAD) from Microorganism 1.2% w/w
Potassium ferricyanide 2.2% w/w
Non-reactive ingredients 96% w/w


References

1. American Diabetes Association, Standards of Medical Care in Diabetes, Diabetes Care, 2016:39 (Supplement 1) S39-S46.
2. Atkins SH, Dasmahapatra A, Jaker MA, Chorost MI, Reddy S: Fingertstick glucose determination in shock. Ann Int Med 114: 1020-1024. 1991
3. Sandler M, Low-Beer T: Misleading capillary glucose measurements. Practical Diabetes 7: 210. 1990
4. Wickham NWR, Achar KN, Cove DH: Unreliability of capillary blood glucose in peripheral vascular disease. Practical Diabetes 3: 100. 1986.

For in vitro diagnostic use

Icon description

- Authorised Representative in the European Community
- Batch code
- Catalogue number
- Consult instructions for use
- Do not reuse
- Do not use if package is damaged
- In vitro diagnostic medical device
- Keep away from sunlight
- Keep dry
- Manufacturer
- Sufficient for
- Temperature limitation
- Use by
- Waste sorting icon
- Imported and distributed by

Spengler SAS
30 rue Jean de Guiramand
13290 Aix-en-Provence - FRANCE

Customer technical support centre:
+33(0)4 42 90 31 31 (press 3)

Tyson Bioresearch, Inc.
5F., No. 16,18,20,22, Kedong 3rd Rd.,
ZhuNan Township, Miaoli County 35053,
Taiwan

Medical Device Safety Service GmbH (MDSS)
Schiffgraben 41, D-30175 Hannover,
Germany

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- American Diabetes Association Clinical Practice Recommendations, Diabetes Care, 2008;31(supplement 1) S12-S54.
- Atkins SH, Dasmahapatra A, Jaker MA, Chorost MI, Reddy S: Fingertstick glucose determination in shock. Ann Int Med 114: 1020-1024. 1991
- Sandler M, Low-Beer T: Misleading capillary glucose measurements. Practical Diabetes 7: 210. 1990
- Wickham NWR, Achar KN, Cove DH: Unreliability of capillary blood glucose in peripheral vascular disease. Practical Diabetes 3: 100. 1986.

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