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JT100 Control Solution Package Insert
Ref : GLU042

INTENDED USE

The Control Solution is a desired aqueous glucose solution with small amounts of non-reactive ingredients and food coloring for verifying the accuracy of the entire JT100 Blood Glucose Monitoring System. The control solutions are used to check the performance of the JT100 meter , JT100 test strips, and your testing technique. The JT100 Blood Glucose Monitoring System is intended for use in the home self-testing and in professional settings to monitor blood glucose levels.

WHEN TO PERFORM A CONTROL SOLUTION TEST

1. You get your meter for the first time before testing your blood test.
2. You open a new vial of test strips.
3. You want to check the meter and strips.
4. You suspect that the meter or test strips are not working properly.
5. Your blood glucose test results are not consistent with how you feel.
6. You think your test result does not accurate.
7. You dropped your meter.
8. Whenever your healthcare professional recommends it.

CONTROL SOLUTION CONTENTS

Each Control Solution is provided as a buffered, stabilized, aqueous solution of desired D(+)-glucose containing a preservative and other non-reactive ingredients. The solutions are formulated to react equivalently to a whole blood sample containing approximate 50, 130 or 350 mg/dL glucose concentrations. Actual lot-specific ranges are printed on the Test Strip bottle label.

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CAUTION: Control Solution is for *in vitro* use only. It is not to be swallowed or injected. The use of the Control Solution is recommended to verify proper system performance and should not be considered as a substitution of actual capillary whole blood testing for blood glucose level.

MATERIALS REQUIRED

1. Glucose Meter
2. Test Strip
3. Control Solution

Note: Please make sure you are using the same model of Test Strip and Glucose Meter.

TEST PROCEDURE

Note: Make sure the solution is at room temperature 4°C to 40°C (39°F to 104°F) prior to testing.

PREPARATION

Make sure that you have set the Glucose Meter in the Control Solution test "[C]" mode, so the result will be stored in the memory and not be incorporated into the average. Please read the User's Manual of your JT100 Blood Glucose Monitoring System.

1. For JT100 Blood Glucose Test Strip, insert a Test Strip, Contact Bar's end first, into the test port. The Meter will automatically turn on. A beep will sound and the system check symbol "[C]" will appear followed by a flashing "" symbol.

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2. After the "" symbol appears, press the M button for 2 sec to switch control solution mode. The "" symbol indicates that the Meter will mark your current test as a Control Solution test. Now you are ready to apply the Control Solution.
3. Discard the first drop of Control Solution and squeeze a small drop of Control Solution on a clean nonabsorbent surface (such as a clean piece of wax paper). Do not apply Control Solution to the Test Strip directly from the bottle.
4. Bring the tip of the Test Strip to lightly touch the drop of Control Solution. Control Solution is automatically pulled into the strip through the tip. Hold until the meter beeps. The meter will now being counting down from 5 to 1 and the Control Solution test result will appear.
5. Compare the results with the Control Solution range printed on the Test Strip vial. The results should fall within the specified range listed on the Test Strip vial.

CAUTION:

- Use JT100 Control Solution, JT100 Blood Glucose Test Strip and JT100 Blood Glucose Meter
- DO NOT apply Control Solution to the Test Strip directly from the bottle, such operation could be leaded to contaminate and decrease the concentration of Control Solution. You MUST place a drop of Control Solution on a clean nonabsorbent surface first.
- Only for use outside the body (For *in vitro* diagnostic use).

COMPARING CONTROL SOLUTION RESULTS

If the test result falls outside the specified range printed on the Test Strip vial, repeat the test.

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RESULTS FALLING OUTSIDE THE SPECIFIED RANGE MAY BE CAUSED BY:

- Error in performing the test.
 - Control solution temperature is lower than 4°C(39°F) or higher than 40°C (104°F).
 - Expired or contaminated Control Solution.
 - Expired or contaminated Test Strips.
 - Improper coding of the Meter.
 - Meter malfunction.
- If the problem persists, the system may not be working properly. **DO NOT** continue to use the system to test your blood glucose level. If you continue to get results that fall outside the specified range, please contact Spengler Customer Technical Support Center : www.sav.gsh-med.fr

STORAGE AND STABILITY

For best results, store Control Solution at room temperature 4°~30°C(39-86°F); do not refrigerate. Refer to the expiration date for stability. The solution should NOT be used after the expiration date or 6 months after opening.

ICON DESCRIPTION

	Authorised Representative in the European Community
	Batch code
	Catalogue number
	Consult instructions for use
	Control
	Do not use if package is damaged
	<i>In vitro</i> diagnostic medical device
	Keep away from sunlight
	Manufacturer
	Temperature limitation
	Use by
	Waste sorting icon
	Imported and distributed by

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**SPENGLER SAS** - 30 rue Jean de Guiramand
**Tyson Bioresearch, Inc.**
**Medical Device Safety Service GmbH (MDSS)**