

ICE FIX TECHNCIAL SPECIFICATION

Definition

- Ice bag.

General features

- Ice bag in tissue comfort double layer
- Hermetic closed tap with safety sheath



Applications

- Unscrew the tap
- Fill in the bag with ice
- Screw the tap and place the bag where necessary



Technical specification

	DIMENSION	
	28 CM	
RAW MATERIAL	Plastic, non-toxic medical level material, biocompatible	
Material used (container):	Terylene	
Material used (tap):	Polythylene	
Nominal Volume (container):	500 ml	
Temperature tollerate (container):	Min -15	
Temperature tollerate (tap):	Min -15	
Dimension container:	Ø 80 x 280 mm	
Dimension tap:	Ø 90 x 20 mm	

Sterile	No
Latex in the product	No
Latex in the package	No
Storage information	Store in a cold dry place away from heat sources

Sterilisation

- Product non sterilisable.

Disposal

- In compliance with the normative in force.

Controlli qualitativi

- Each batch is controlled according to the internal procedures. Controls are made on incoming raw material and during the manufacturing process. Before proceeding with the packaging are effected some aesthetical and visual and dimensional tests.



ICE FIX TECHNICAL SPECIFICATION

Package

Labels and/or packages have the following:

Name of sell
Item barcode

Packaging material

- Box in carton Pic branded.

Supplier

- Artsana S.p.A. Via Saldarini Catelli 1, 22070 Grandate (Como) Italy.

Assortment

Code	Box	Unit	Package
00 012302 100 000	1 pz	6 pz	24 pz (4 units)

The Technical Department

ARTSANA S.p.A.
Dr. MAURO CASSANI
Medical Devices Department
Technical Director

Artsana S.p.A. considers the information in this document confidential and non-disclosable. Artsana reserves all rights to change this document based on technical developments and amendments to relative product standards.

