



## **TECHNICAL DATA SHEET** **ICE BAGS IN PE**

**NAME OF THE PRODUCT:** EASY ICE – DISPO ICE – MINI DISPO ICE – IMPLANT ICE

**DESCRIPTION OF THE PRODUCT:** INSTANT ICE BAGS

**DIMENSIONS:** Cm 14x18, 14x24, tolerance +/-5%

**CLASSIFICATION:** Medical Device Class IIA

Classification according to the Annex IX paragraph III rule 9 of the European Directive 93/42/EEC (and successive amendments included Directive 2007/47/EEC)

### **TECHNICAL FEATURES:**

**External bag:** PE with three heat seals so to assure thermo soldering, impermeableness and resistance. PE: 2 layers.

**Internal water bag** made of special 2-layer material so to break with a slight pressure of the fingers, shockproof. Covered by patent.

**Content:** urea, ammonium chloride and water.

**MANUFACTURER:** Dispotech srl Via M. del Grosso 8-12  
23022 Chiavenna (SO)  
Tel. 0343.36711  
Fax 0343.36567  
e-mail: [info@dispotech.it](mailto:info@dispotech.it)

**TEMPERATURE:** Thanks to the mixture, the temperature reaches -4°C arriving to 0°C in 20 minutes

**LOT:** Mentioned on the back of the bag

**EXPIRATION DATE:** 5 years from production date. It is mentioned on the back of the bag

**SYMBOLS:** Symbols according to European rule UNI CEI EN ISO 15223



**STORAGE:** Keep at ambient temperature, in a dry place, away from heat sources, flames or sparks.

**PACKAGING:** boxes of 25pcs

DATE: 23/11/2018  
Rev. 06  
MADE BY RGQ  
APPROVED BY DG

### **Dispotech srl**

Sede operativa: 23020 GORDONA (SO) Via al piano, 29 - Tel. +39 0343 36711 - Fax. +39 0343 36567 - [www.dispotech.com](http://www.dispotech.com)

Sede legale: 23022 CHIAVENNA (SO) Via M. Del Grosso, 8/12 - e-mail: [info@dispotech.it](mailto:info@dispotech.it) - PEC: [dispotech@legalmail.it](mailto:dispotech@legalmail.it)

Cap Soc. € 1.500.000,00 i.v. - P. IVA 00672170149 - R.E.A. 47213 C.C.I.A.A. di SO - Uff. Reg. Imp. SO 00672170149

