

Declaration of Conformity

In accordance with Annex II Clause 3 of the MDD 93/42/EEC

Manufacturer: UAB VITACON LT
Address: Naujoji g.12-525.
Alitaus 62401
Lithuania

PRODUCT:

Type/Model: VitaScan PD Ultrasound Bladder Scanner
Description: Unique non-invasive system for automatic bladder volume measurement.
Classification: IIa

VitaScan PD Ultrasound Bladder Monitor

Name as per Device Label	Identifier (Item Code)	Brief Description of Item
VitaScan PD Bladder Scanner	100550	The VitaScan PD is a B-mode ultrasonic instrument, which projects ultrasound energy through the lower abdomen of the non-pregnant patient to obtain an image of the bladder, which is used to determine bladder volume non-invasively, on an intermittent basis. Including: VitaScan PD console and USB probe, Probe Sleeve, VitaScan PD software version 1.0.3, Digital documentation, and Calibration certificate.
VitaScan PD Bladder Scanner - Medical grade system with VitaScan Mobile Cart	100550 MP	Including: VitaScan PD Bladder Scanner - USB probe, VitaScan console with preinstalled VitaScan PD software version 1.0.3, VESA mount Bracket with adhesive tape, VitaScan Mobile Cart with adjustable height, VitaScan Probe Sleeve, Probe Brace and Basket with handle, Digital documentation.
VitaScan Mobile Cart	100762	VitaScan Mobile Cart with adjustable height, Probe Brace and Basket with handle
VitaScan Calibration Test Tool	100651	Test Tool to perform VitaScan PD on-site calibration
VitaScan Verification Test Tool	100752	Test Tool to perform VitaScan PD on-site verification
VitaScan PD console	100551	Console of VitaScan PD Bladder scanner, excluding the Probe
VitaScan replacement probe	100664	Replacement probe for VitaScan Bladder scanner

We herewith declare that above mentioned product meet the transposition into national law of the provisions of council directive MDD 93/42/EEC

UAB VITACON LT

Date: 03rd of November, 2017

Signature: Svein Wognild
Svein Wognild

VITACON.