

EAR & FOREHEAD INFRARED THERMOMETER

CT5D1



Feature

- Applicable for forehead and ear temperature
- Special ergonomic design to facilitate operation
- Memory function recalls 20 previous readings at most
- No danger of breaking the glass or swallowing mercury
- Fever warning when the reading exceeds 37.5 °C / 99.5 °F
- Large digital LCD display with backlight
- CE mark and FDA approved

Specifications

Measurement range	Human body model: 32-43 °C / 89.6-109.4 °F The object model: 0-60 °C / 32-140 °F		
Auto shut-off function	Within 60s	Minimum scale	0.1 °C / °F
Display resolution	0.1 °C / 0.1 °F	Accuracy	±0.2 °C / ±0.4 °F
Dimensions	146(L)*43.8(W)*39.9(H) mm	Weight	70.7g
Battery life	More than 3000 takes	Battery	2 AAA alkaline batteries
Standard accessories	1 user manual; 2 AAA batteries	Applicable user	Adults and children