



PESE-BEBE « CARTOON DESIGN » + SACOCHE



CARACTERISTIQUES :

- Portée maximale : 20Kg
- Graduation : 10g. de 0 à 20kg
- Grand afficheur LCD rétro éclairé 65 X 30 mm
- Equipé de 4 capteurs « haute précision »
- Nacelle extra large permettant le pesage assis ou allongé
- Possibilité de mélodie musicale pendant la pesée
- Mémorisation du dernier poids
- Boutons tare, unité de mesure, musique
- Tare automatique
- Alimentation : 4 piles LR6 fournies
- Indicateur d'usure des piles
- Accessoires fournis : Sacoche, mètre ruban, alèse
- Poids de l'appareil : 1kg600 Dimension : 62 x 32 x 7 cm
- Coloris Blanc : Référence P25 BLANC, Code Article : 3 700445 50024 6
- Coloris Bleu : Référence P25 BLEU, Code Article : 3 700445 50026 0



Carton de 2 pcs
Carton de 2 pcs

The image shows the product packaging for the LBS Medical baby scale. It features a central image of the white scale and a list of accessories: a blue carrying bag, a measuring tape, a scale mat, and a manual. The packaging is multilingual, with text in French, English, Dutch, and Spanish. The French text reads: 'PÈSE BÉBÉ ÉLECTRONIQUE (Mode d'emploi détaillé à l'intérieur)'. The English text reads: 'ELECTRONIC BABY SCALE (Detailed manual inside)'. The Dutch text reads: 'BABYWEEGSCHAAL (Uitgebreide handleiding bij de binnenzijde)'. The Spanish text reads: 'PESA BÉBÉ ELECTRÓNICA (Manual detallado en el interior)'. There are also small icons for a baby, a 20kg weight limit, and a battery symbol.

LBS MEDICAL

LBS 101 CHEMIN DE RONDE 78290 CROISSY SUR SEINE TEL : 01 30 15 06 96 FAX : 01 39 76 45 52

E-mail : labo.lbs@free.fr www.lbs-medical.com