

optimized

seca 336

Electronic baby scale



- Handy and light for mobile medical use.
- Flat cradle.
- Highly legible, illuminated LCD display.
- Optional digital measuring rod seca 234 or analog measuring rod seca 232 n.

seca 336:

High quality and convenient operation
for fast, accurate results.

► Modern technology is helpful only if it works in practice. The electronic baby scale seca 336 has all the features needed to ease the everyday workload in medical practices and hospitals. Simple to use, mobile and durable, the scale satisfies all the requirements of an ideal scale in pediatrics.



Lightweight and handy for
mobile medical use.

Weighing just 3.7 kg, the scale is very lightweight in its class. The handle integrated on the side makes transport easy. Another practical aspect is its flat and space-saving shape, which fits into the drawer of a changing table. For mobile use the scale is safely transported in the optionally available carrying case seca 428.



With or without measuring rod, the seca 336 is suitable for transport, thanks to the handle and its own low weight.

Convenient operation and
easy-to-read LCD display.

The spacious, flat cradle gives baby comfort and security. With the little one comfortably positioned on the cradle, the length measurement can be carried out quickly and precisely. Larger babies or small children can sit while being weighed. The measurement is shown on the brightly lit LCD display, where it remains for up to two minutes.



Easy-to-read, illuminated LCD display.

A faster measuring process is
more pleasant for the baby.

The seca 336 is equipped with a special damping system that delivers exact measurements immediately even when baby is restless. Combined with the optional measuring rod, either the digital seca 234 or analog seca 232 n, the scale lets you measure weight and length in one step.



The digital measuring rod with a clear display simplifies the measuring process.

Technische Daten

- Capacity: 20 kg / 44 lbs
- Graduation: 5 g < 10 kg > 10 g / 0.2 oz < 22 lbs > 0.5 oz
- Dimensions (WxHxD): 650 x 110 x 308 mm / 25.6 x 4.3 x 12.1"
- Weight: 3.7 kg / 7.7 lbs
- Functions: TARE, BMIF, HOLD, Auto-HOLD, auto weighing range switch-over, damping, automatic switch-off, Confirm
- Optional: mobile carts seca 403 and seca 402, digital measuring rod seca 234, analog measuring rod seca 232 n, power adapter, carrying case seca 428