



STE 292-H L

Dual Motor Lift Chair

with Headrest and Lumbar System

COMFORT & SMART FEATURES

Dual Motor Lift & Layflat Mechanism

Steel Seat Box

Footrest Extension

Power Headrest and Lumbar System

**Chaise Pad Design Fully Supports Legs
Eliminating Pressure Points**

**Handwand with Independent
Backrest & Footrest Control**

All Close and Lift Button

USB Charging Port

SPECIFICATIONS

Depth / Width / Height: 40.5" / 28.8" / 38"

Seat Depth / Width / Height: 18.5" / 19.5" / 18"

Fully Reclined Length: 70"

Wall Clearance Required: 16"

Foam Density: 2.2

Weight Capacity: 300 lbs

Loading QTY / 40' HC: 144

Cubes: 15.7

Box Dimensions: 28.3" x 30.3" x 31.5"

Net Weight: 124.5 lbs

Gross Weight: 142.2 lbs

Cover # 1		Cover # 2	
Color		Color	
Number		Number	
Price		Price	