


Adjustable walking cane Instructions

Thank you for purchasing our Adjustable Walking Stick, Our adjustable walking cane will assist you walking easiler. If you want to adjust the height, please press the buttons then you will adjust your suitable height.



- HOW TO USE IT -



Going Upstairs Support



Going Downstairs Support




Standing Up Support

Canne de marche ajustable

Mode d'emploi

Merci d'avoir acheté notre de canne de marche ajustable, elle vous aidera à vous déplacer plus facilement.

Pour ajuster la hauteur de la canne, appuyez sur le bouton sous le pommeau avec une main tout en tirant le corps avec l'autre jusqu'à la hauteur désirée.



Exemples d'utilisation



MONTER LES ESCALIERS



DESCENDRE LES ESCALIERS



AIDE AU LEVER

